


4ME 16:0 PA sodium

Chemical Properties

CAS No. :	474967-75-4
Formula:	C43H84NaO8P
Molecular Weight:	783.09
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	4ME 16:0 PA sodium is a bis(monoacylglycero)phosphate lipid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.277 mL	6.385 mL	12.7699 mL
5 mM	0.2554 mL	1.277 mL	2.554 mL
10 mM	0.1277 mL	0.6385 mL	1.277 mL
50 mM	0.0255 mL	0.1277 mL	0.2554 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481