

20:4-16:0 PC - IsoPure

## Chemical Properties


CAS No. : 118139-93-8

Formula: C44H80NO8P

Molecular Weight: 782.08

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	20:4-16:0 PC-IsoPure is a mixed acylglycerophospholipid.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2786 mL	6.3932 mL	12.7864 mL
5 mM	0.2557 mL	1.2786 mL	2.5573 mL
10 mM	0.1279 mL	0.6393 mL	1.2786 mL
50 mM	0.0256 mL	0.1279 mL	0.2557 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481