


Taurodehydrocholic acid sodium

Chemical Properties

CAS No. :	57011-24-2
Formula:	C ₂₆ H ₃₈ NNaO ₇ S
Molecular Weight:	531.64
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Taurodehydrocholic acid sodium acts as an activator of biliary cholesterol secretion. It notably increases the expression of Abcg5 while decreasing the expression of abc8. This compound is applicable in research related to cholesterol metabolism.
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.881 mL	9.4049 mL	18.8097 mL
5 mM	0.3762 mL	1.881 mL	3.7619 mL
10 mM	0.1881 mL	0.9405 mL	1.881 mL
50 mM	0.0376 mL	0.1881 mL	0.3762 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481