

C25 Ceramide (d18:1/25:0)

Chemical Properties


CAS No. : 67492-18-6

Formula: C₄₃H₈₅NO₃

Molecular Weight: 664.14

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	C25 Ceramide (d18:1/25:0) is a type of ceramide.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5057 mL	7.5285 mL	15.0571 mL
5 mM	0.3011 mL	1.5057 mL	3.0114 mL
10 mM	0.1506 mL	0.7529 mL	1.5057 mL
50 mM	0.0301 mL	0.1506 mL	0.3011 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481