


Methyl N-boc-3-chloro-L-alaninate

Chemical Properties

CAS No. : 651035-84-6
 Formula: C₉H₁₆ClNO₄
 Molecular Weight: 237.68
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Methyl N-boc-3-chloro-L-alaninate is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2073 mL	21.0367 mL	42.0734 mL
5 mM	0.8415 mL	4.2073 mL	8.4147 mL
10 mM	0.4207 mL	2.1037 mL	4.2073 mL
50 mM	0.0841 mL	0.4207 mL	0.8415 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481