


Chenodeoxycholic acid disulfate sodium

Chemical Properties

CAS No. :	71387-96-7
Formula:	C ₂₄ H ₃₇ Na ₃ O ₁₀ S ₂
Molecular Weight:	618.64
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Chenodeoxycholic acid disulfate sodium is a derivative of Chenodeoxycholic acid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6164 mL	8.0822 mL	16.1645 mL
5 mM	0.3233 mL	1.6164 mL	3.2329 mL
10 mM	0.1616 mL	0.8082 mL	1.6164 mL
50 mM	0.0323 mL	0.1616 mL	0.3233 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481