


trans-3-Hexenoyl-CoA

Chemical Properties

CAS No. :	157893-61-3
Formula:	C ₂₇ H ₄₄ N ₇ O ₁₇ P ₃ S
Molecular Weight:	863.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	trans-3-Hexenoyl-CoA (3-trans-Hexenoylcoenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1579 mL	5.7893 mL	11.5786 mL
5 mM	0.2316 mL	1.1579 mL	2.3157 mL
10 mM	0.1158 mL	0.5789 mL	1.1579 mL
50 mM	0.0232 mL	0.1158 mL	0.2316 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481