


18:1(n9)-16:0 PG-IsoPure sodium

Chemical Properties

CAS No. :	2389048-55-7
Formula:	C40H76NaO10P
Molecular Weight:	770.99
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	18:1(n9)-16:0 PG-IsoPure sodium is a mixed acyl phospholipid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.297 mL	6.4852 mL	12.9703 mL
5 mM	0.2594 mL	1.297 mL	2.5941 mL
10 mM	0.1297 mL	0.6485 mL	1.297 mL
50 mM	0.0259 mL	0.1297 mL	0.2594 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481