

## trans-2-Methyl-5-isopropylhexa-2,5-dienoyl-CoA

## Chemical Properties


CAS No. : 1245945-94-1

Formula: C<sub>31</sub>H<sub>50</sub>N<sub>7</sub>O<sub>17</sub>P<sub>3</sub>S

Molecular Weight: 917.75

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	trans-2-Methyl-5-isopropylhexa-2,5-dienoyl-CoA (trans-2-Methyl-5-isopropylhexa-2,5-dienoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0896 mL	5.4481 mL	10.8962 mL
5 mM	0.2179 mL	1.0896 mL	2.1792 mL
10 mM	0.109 mL	0.5448 mL	1.0896 mL
50 mM	0.0218 mL	0.109 mL	0.2179 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481