

## Photoclick C6 ceramide

## Chemical Properties


CAS No. : 2389048-46-6

Formula: C<sub>24</sub>H<sub>41</sub>N<sub>3</sub>O<sub>3</sub>

Molecular Weight: 419.60

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Photoclick C6 ceramide is a clickable sphingolipid.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3832 mL	11.9161 mL	23.8322 mL
5 mM	0.4766 mL	2.3832 mL	4.7664 mL
10 mM	0.2383 mL	1.1916 mL	2.3832 mL
50 mM	0.0477 mL	0.2383 mL	0.4766 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481