


## 18:1 PE-Benzylguanine ammonium

## Chemical Properties

CAS No. :	2315262-37-2
Formula:	C <sub>55</sub> H <sub>93</sub> N <sub>8</sub> O <sub>10</sub> P
Molecular Weight:	1057.35
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	18:1 PE-Benzylguanine ammonium is a lipid that reacts with the SNAP tag.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9458 mL	4.7288 mL	9.4576 mL
5 mM	0.1892 mL	0.9458 mL	1.8915 mL
10 mM	0.0946 mL	0.4729 mL	0.9458 mL
50 mM	0.0189 mL	0.0946 mL	0.1892 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481