

PFVYLI

Chemical Properties

CAS No. : 1252672-28-8


Formula: C40H58N6O8

Molecular Weight: 750.92

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PFVYLI is a hydrophobic penetrating peptide used to modify stealth sustained-release liposomes loaded with Doxorubicin (PFV-ssls-dox) for breast cancer research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3317 mL	6.6585 mL	13.317 mL
5 mM	0.2663 mL	1.3317 mL	2.6634 mL
10 mM	0.1332 mL	0.6658 mL	1.3317 mL
50 mM	0.0266 mL	0.1332 mL	0.2663 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481