

28:0 Coenzyme A ammonium

Chemical Properties


CAS No. : 2260795-80-8

Formula: C₄₉H₉₉N₁₀O₁₇P₃S

Molecular Weight: 1225.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	The compound 28:0 Coenzyme A ammonium is a saturated fatty acid acyl-CoA.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8161 mL	4.0805 mL	8.1609 mL
5 mM	0.1632 mL	0.8161 mL	1.6322 mL
10 mM	0.0816 mL	0.408 mL	0.8161 mL
50 mM	0.0163 mL	0.0816 mL	0.1632 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481