


S-(Hydrogen3-ethyl-3-hydroxypentanedioate)-CoA

Chemical Properties

CAS No. :	69925-08-2
Formula:	C ₂₈ H ₄₆ N ₇ O ₂₀ P ₃ S
Molecular Weight:	925.69
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	S-(Hydrogen3-ethyl-3-hydroxypentanedioate)-CoA (Coenzyme A, S-(hydrogen3-ethyl-3-hydroxypentanedioate)) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0803 mL	5.4014 mL	10.8028 mL
5 mM	0.2161 mL	1.0803 mL	2.1606 mL
10 mM	0.108 mL	0.5401 mL	1.0803 mL
50 mM	0.0216 mL	0.108 mL	0.2161 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481