


16:0 Hexynoyl PE ammonium

Chemical Properties

CAS No. :	2260670-40-2
Formula:	C43H83N2O9P
Molecular Weight:	803.10
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	16:0 Hexynoyl PE ammonium is a clickable phospholipid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2452 mL	6.2259 mL	12.4517 mL
5 mM	0.249 mL	1.2452 mL	2.4903 mL
10 mM	0.1245 mL	0.6226 mL	1.2452 mL
50 mM	0.0249 mL	0.1245 mL	0.249 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481