


## trans-2,5-Diethylpiperazine

## Chemical Properties

CAS No. :	6189-24-8
Formula:	C <sub>8</sub> H <sub>18</sub> N <sub>2</sub>
Molecular Weight:	142.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	trans-2,5-Diethylpiperazine is a pharmaceutical intermediate used in the synthesis of various active compounds.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.0304 mL	35.1519 mL	70.3037 mL
5 mM	1.4061 mL	7.0304 mL	14.0607 mL
10 mM	0.703 mL	3.5152 mL	7.0304 mL
50 mM	0.1406 mL	0.703 mL	1.4061 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481