


## GT3 Sphingosine (d18:1)

## Chemical Properties

CAS No. :	736131-61-6
Formula:	C63H108N4O36
Molecular Weight:	1497.54
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	GT3 Sphingosine (d18:1) is a type of ganglioside.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6678 mL	3.3388 mL	6.6776 mL
5 mM	0.1336 mL	0.6678 mL	1.3355 mL
10 mM	0.0668 mL	0.3339 mL	0.6678 mL
50 mM	0.0134 mL	0.0668 mL	0.1336 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481