


## Butyl diphenyl phosphate

## Chemical Properties

CAS No. :	2752-95-6
Formula:	C <sub>16</sub> H <sub>19</sub> O <sub>4</sub> P
Molecular Weight:	306.29
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Butyl diphenyl phosphate is primarily utilized as a flame retardant and plastic additive.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2649 mL	16.3244 mL	32.6488 mL
5 mM	0.653 mL	3.2649 mL	6.5298 mL
10 mM	0.3265 mL	1.6324 mL	3.2649 mL
50 mM	0.0653 mL	0.3265 mL	0.653 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481