


## Sulfoacetyl-CoA

## Chemical Properties

CAS No. : 1260155-48-3  
 Formula: C<sub>23</sub>H<sub>38</sub>N<sub>7</sub>O<sub>20</sub>P<sub>3</sub>S<sub>2</sub>  
 Molecular Weight: 889.63  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sulfoacetyl-CoA (Sulfoacetyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1241 mL	5.6203 mL	11.2406 mL
5 mM	0.2248 mL	1.1241 mL	2.2481 mL
10 mM	0.1124 mL	0.562 mL	1.1241 mL
50 mM	0.0225 mL	0.1124 mL	0.2248 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481