

## N-Acetylalanine-CoA

## Chemical Properties


CAS No. : 1048681-16-8

Formula: C<sub>26</sub>H<sub>43</sub>N<sub>8</sub>O<sub>18</sub>P<sub>3</sub>S

Molecular Weight: 880.65

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-Acetylalanine-CoA (N-Acetylalanine-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1355 mL	5.6776 mL	11.3552 mL
5 mM	0.2271 mL	1.1355 mL	2.271 mL
10 mM	0.1136 mL	0.5678 mL	1.1355 mL
50 mM	0.0227 mL	0.1136 mL	0.2271 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481