

CP-LC-1421

## Chemical Properties


CAS No. : 3091879-18-1

Formula: C<sub>57</sub>H<sub>109</sub>N<sub>3</sub>O<sub>8</sub>S<sub>3</sub>

Molecular Weight: 1060.69

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CP-LC-1421 is a certest lipid.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9428 mL	4.7139 mL	9.4278 mL
5 mM	0.1886 mL	0.9428 mL	1.8856 mL
10 mM	0.0943 mL	0.4714 mL	0.9428 mL
50 mM	0.0189 mL	0.0943 mL	0.1886 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481