


## Anhydrolutein II

## Chemical Properties

CAS No. :	92760-19-5
Formula:	C <sub>40</sub> H <sub>54</sub> O
Molecular Weight:	550.86
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Anhydrolutein II is a type of carotenoid and is the dehydrated derivative of lutein.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8153 mL	9.0767 mL	18.1534 mL
5 mM	0.3631 mL	1.8153 mL	3.6307 mL
10 mM	0.1815 mL	0.9077 mL	1.8153 mL
50 mM	0.0363 mL	0.1815 mL	0.3631 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481