

C11 BODIPY 505/515 MG

Chemical Properties

CAS No. : 2260795-64-8


Formula: C27H41BF2N2O4

Molecular Weight: 506.43

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	C11 BODIPY 505/515 MG is a glycerolipid tagged with a fluorescent marker (Ex/Em = 495/503 nm).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9746 mL	9.873 mL	19.7461 mL
5 mM	0.3949 mL	1.9746 mL	3.9492 mL
10 mM	0.1975 mL	0.9873 mL	1.9746 mL
50 mM	0.0395 mL	0.1975 mL	0.3949 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481