


BODIPY 505/515 TMR PI(3,4,5)P3 ammonium

Chemical Properties

CAS No. :	3077081-40-1
Formula:	C54H97BF2N7O24P4
Molecular Weight:	1401.08
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	BODIPY 505/515 TMR PI(3,4,5)P3 ammonium is a fluorescently labeled phospholipid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7137 mL	3.5687 mL	7.1374 mL
5 mM	0.1427 mL	0.7137 mL	1.4275 mL
10 mM	0.0714 mL	0.3569 mL	0.7137 mL
50 mM	0.0143 mL	0.0714 mL	0.1427 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481