

## BDP TMR PI(3)P diammonium

## Chemical Properties

CAS No. : 3055900-46-1


Formula: C54H88BF2N5O18P2

Molecular Weight: 1206.06

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BDP TMR PI(3)P diammonium is a phospholipid with a fluorescent label.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8291 mL	4.1457 mL	8.2915 mL
5 mM	0.1658 mL	0.8291 mL	1.6583 mL
10 mM	0.0829 mL	0.4146 mL	0.8291 mL
50 mM	0.0166 mL	0.0829 mL	0.1658 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481