


Hexyl 4-hydroxybenzoate

Chemical Properties

CAS No. :	1083-27-8
Formula:	C ₁₃ H ₁₈ O ₃
Molecular Weight:	222.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hexyl 4-hydroxybenzoate (Hexylparaben) is an agonist of the retinoic acid-related orphan receptor ROR γ (t) with an EC ₅₀ of 144 nM. It is applicable in the study of autoimmune diseases.
Targets(IC ₅₀)	ROR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4988 mL	22.4942 mL	44.9883 mL
5 mM	0.8998 mL	4.4988 mL	8.9977 mL
10 mM	0.4499 mL	2.2494 mL	4.4988 mL
50 mM	0.090 mL	0.4499 mL	0.8998 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481