


SNAP-fusion protein ligand-1

Chemical Properties

CAS No. :	913817-45-5
Formula:	C ₁₂ H ₁₃ ClN ₄ O
Molecular Weight:	264.71
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	SNAP-fusion protein ligand-1 (Compound 1i) serves as a ligand for a target protein in PROTAC (Ligand for Target Protein for PROTAC) applications. It is used in the synthesis of PROTAC VHL-SNAP2-5C.
Targets(IC50)	Others,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7777 mL	18.8886 mL	37.7772 mL
5 mM	0.7555 mL	3.7777 mL	7.5554 mL
10 mM	0.3778 mL	1.8889 mL	3.7777 mL
50 mM	0.0756 mL	0.3778 mL	0.7555 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481