

50308,12

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	50308,12 is an ionizable cationic lipidoid utilized in forming lipid nanoparticles (LNPs) for mRNA delivery in living organisms. LNPs constructed with 50308,12 that encapsulate an mRNA reporter are known to target the spleen in mice.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481