

8-Arm-PEG-FA (MW 600)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	8-Arm-PEG-FA (MW 600) [8-Arm-PEG-folate (MW 600)] is a PEG derivative where each terminal of the PEG chains is linked to a Folic acid molecule. Folic acid has a high affinity for folate receptors and can be utilized for targeting cell membrane receptors in drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481