

Alkyne-PEG-COOH (MW 40000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Alkyne-PEG-COOH (MW 40000) is a heterobifunctional PEG cross-linker that facilitates the attachment of functional PEG to biomolecules, particles, and other material surfaces. The carboxyl group can react with amine groups to form a stable amide bond.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481