

DSPE-PEG3400-YIGSR

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG3400-YIGSR is a PEG compound composed of DSPE and the biomimetic peptide YIGSR. YIGSR interacts with the 67 kDa laminin binding protein (LBP), promoting adhesion and proliferation of various cell types, including endothelial cells, fibroblasts, and smooth muscle cells. DSPE-PEG3400-YIGSR can be utilized in drug delivery applications.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481