

15-Deoxy- Δ 12,14-prostaglandin J2 serinol amide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	15-Deoxy- Δ 12,14-prostaglandin J2 serinol amide (15-Deoxy- Δ 12,14-PGJ2 serinol amide) is a derivative of 15-deoxy- Δ 12,14-PGJ2-2-glycerol ester.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481