

## Urocortin II, human acetate

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Urocortin II, human acetate, is a selective endogenous peptide agonist of type 2 corticotropin-releasing factor (CRF2) receptors, used to study the role of CRF2 receptors in feeding behavior.
Targets(IC50)	CRFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481