

4-Arm-PEG-FA (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | 4-Arm-PEG-FA (MW 2000) [4-Arm-PEG-folate (MW 2000)] is a PEG derivative with each PEG chain terminating in a folic acid molecule. Folic acid has a high affinity for folate receptors, making it useful for targeting cell membrane receptors in drug delivery. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481