

## GlyCOUPER protease

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	GlyCOUPER protease is suitable for digesting fusion protein flexible linkers composed of glycine or glycine and serine residues, such as (Gly4Ser) <sub>n</sub> , GlyxSery (GS), and polyglycine (G) linkers.
Targets(IC50)	Serine Protease

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481