

## DOTA-Bombesin (1-14) TFA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DOTA-Bombesin (1-14) TFA is composed of the chelating agent DOTA and Bombesin. It is utilized in cancer research by binding with radioactive elements. DOTA-Bombesin (1-14) TFA can be applied in the synthesis and study of radiopharmaceutical conjugates
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481