

DSPE-PEG3400-GE11

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | DSPE-PEG3400-GE11 is a PEG compound composed of DSPE and an EGFR-targeting peptide (GE11). GE11 is applicable for cancer cells with EGFR overexpression. DSPE-PEG3400-GE11 can be used for drug delivery. |
| Targets(IC50) | EGFR,Liposome |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481