

Boc-NH-PEG-COOH (MW 40000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Boc-NH-PEG-COOH (MW 40000) is a linear heterobifunctional polyethylene glycol (PEG) compound, featuring a Boc-protected amine and a carboxyl group. It serves as a significant crosslinking agent with a PEG chain.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481