

## Biotin-PEG-Biotin(MW 40000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Biotin-PEG-Biotin (MW 40000) is used to crosslink and pegylate by binding with two streptavidin and avidin molecules. The biotin is linked to linear polyethylene glycol (PEG) via a stable amide bond.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481