

DSPE-PEG3400-ANG

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | DSPE-PEG3400-ANG is a PEG compound formed from DSPE and the dual-targeting ligand (Angiopep-2, ANG). ANG demonstrates high binding efficiency to LRP1 and has been employed for targeted delivery to gliomas. DSPE-PEG3400-ANG is suitable for drug delivery applications. |
| Targets(IC50) | Liposome |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481