

BOC-NH-PEG-NH2 (MW 5000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BOC-NH-PEG-NH2 (MW 5000) (Boc-NH-PEG-Amine (MW 5000)) is a linear heterobifunctional PEG compound, featuring one Boc-protected amine and one free amino group. This compound serves as a crucial crosslinker with a PEG chain.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481