

(R,R)-dWIZ-1 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	(R,R)-dWIZ-1 TFA is the (R,R) enantiomer of dWIZ-1. dWIZ-1 functions as a potent molecular glue degrader for WIZ and shows potential in sickle cell disease (SCD) research.
Targets(IC50)	Histone Acetyltransferase,Molecular Glues

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481