

HOOC-PEG-SH (MW 20000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	HOOC-PEG-SH (MW 20000) [HOOC-PEG-Thiol (MW 20000)] is a reactive thiol PEG derivative featuring a terminal carboxyl group. The carboxyl group can react with amines or hydroxyl groups to form either stable amide bonds or unstable ester bonds. The PEG linkage between the thiol and carboxyl groups offers excellent water solubility, flexible joint distance, and enhanced stability.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481