

Biotin-PEG-SH (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Biotin-PEG-SH (MW 3400) (Biotin-PEG-Thiol (MW 3400)) is a biotin-conjugated PEG polymer containing a thiol (-SH) group. Compounds labeled with biotin can bind to avidin or streptavidin for subsequent purification or detection.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481