

mPEG-CM (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	mPEG-CM (MW 1000), also known as mPEG-COOH (MW 1000), is a compound featuring a reactive carboxyl group (-COOH) at the terminal end. This functional group allows it to form stable amide bonds with amino groups or ester bonds with hydroxyl groups. mPEG-CM (MW 1000) is applicable in drug delivery systems.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481