

## 4-Arm PEG-SH (MW 5000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	4-Arm PEG-SH (MW 5000) (4-Arm PEG-thiol (MW 5000)) is a multi-arm PEG derivative with a thiol group at the end of each of its four arms. The reactive free thiol, SH, or mercapto groups selectively react with maleimide and surfaces of transition metals including gold and silver.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481