

DSPE-PEG3400-cRGD

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG3400-cRGD is a PEG compound composed of DSPE and the $\alpha\beta 3$ -targeting peptide (cRGD). The cRGD peptide specifically binds to the $\alpha\beta 3$ integrins present on the surface of various cancer cells and neovascularization cells. DSPE-PEG3400-cRGD can be utilized for drug delivery purposes.
Targets(IC50)	Integrin,Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481