

DSPE-PEG3400-M2pep

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG3400-M2pep is a PEG compound composed of DSPE and the M2 macrophage-targeting peptide (M2pep). M2pep acts as a monovalent or multivalent ligand conjugated with pro-apoptotic peptides to enhance the toxic effect on M2 macrophages. DSPE-PEG3400-M2pep is applicable for drug delivery.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481