

DSPE-PEG5000-PTP

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG5000-PTP is synthesized by covalently linking the amino group of plectin-1 peptide (PTP) to DSPE-PEG-NHS ester through an amide bond. DSPE-PEG1000-PTP can be utilized in the construction of drug delivery carriers.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481