

DSPE-PEG-CHO ammonium (MW 40000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG-CHO ammonium (MW 40000) is a linear heterobifunctional polyethyleneglycolated reagent combining DSPE phospholipids with an aldehyde group. The PEG linker offers excellent hydrophilicity and water solubility. Aldehyde polyethylene glycol (DSPE) is a commonly used reactive phospholipid for conjugating antibodies, peptides, or other ligands to the surface of liposomes and lipid PEG nanoparticles. DSPE-PEG-CHO ammonium (MW 40000) is applicable for drug delivery.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481